


Precision Balances, Nimbus®



0.0001 g Balance

0.001 g Balance

- Vivid, backlit LCD and sealed, color-coded keypad
- Level indicator & adjusting feet ensure proper balance setup
- Security slot for optional Kensington-type lock & cable to help prevent theft
- Removable draft shield on models 0.001g readability
- Robust metal housing protects internal components
- Operating Temperature: 0° to 40°C / 32° to 104°F
- **Units of measure: 0.0001 g Readability Balances:** g, mg, ct, GN, oz, ozt, dwt, custom unit
- **0.001 g Readability Balances:** g, kg, ct, GN, N, lb, oz, ozt, dwt, custom unit
- USB and RS-232 interface
- External calibration allows for verification and adjustment with weights
- Printouts include date and time for data tracking within GLP guidelines
- Selectable digital filtering helps minimize effects of vibration and disturbances
- Zero-tracking feature ensures display returns to zero reading
- 18 VDC 800mA adapter
- AC Adapter included
- 1 YEAR WARRANTY

	Cap. (g)	Read. (g)	Repeatability (std g)	Linearity (g)	Stabilization Time (sec.)	Draft Shield (mm)
M113482	80	0.0001	±0.00015	±0.0002	3	Chamber 165 x 145 x 240
M113483	120	0.0001	±0.00015	±0.0002	3	Chamber 165 x 145 x 240
M113484	160	0.0001	±0.0002	±0.0002	3	Chamber 165 x 145 x 240
M113485	210	0.0001	±0.0002	±0.0002	3	Chamber 165 x 145 x 240
M113486	250	0.0001	±0.0002	±0.0002	3	Chamber 165 x 145 x 240
M113487	220	0.001	±0.002	±0.002	3	Round 180 ø x 110
M113488	420	0.001	±0.002	±0.002	3	Round 180 ø x 110
M113489	620	0.001	±0.002	±0.002	3	Round 180 ø x 110
M113490	820	0.001	±0.002	±0.002	3	Round 180 ø x 110