

Microcentrifuge, Spectrafuge™ 24D



- Precise microprocessor control
- “MultiFlo” airflow pattern keeps samples cool
- Powerful, maintenance-free brushless motor, quiet operation
- Simple knob adjustment maintains precise control; digital display
- Solid aluminum rotor (24 x 1.5/2.0ml included) is easy to remove for cleaning and autoclaving
- Optional StripSpin™ adapter for PCR strips or 16 individual 0.2µl tubes

	M107409	M107410	M107411	M107412	M107413
Color:	Grey	Blue	Teal	Purple	Red
Speed:	500 to 13300 rpm / 16300 x g				
Timer:	1 to 30 mins. or continuous “Quick” Button for momentary operation				
Dims./Wgt.:	9.25 (W) x 11.5 (D) x 8.5 (H) in. / 18 lbs.				



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002