

Mini Centrifuge, My Fuge™ Mini



- Eight Position Microtube rotor 8 x 1.5/2.0ml (comes installed)
- Rotor for 0.2ml PCR tubes & strips (included) is stored in a convenient compartment on bottom of unit
- Near silent operation
- Starts and stops in seconds with open/close of lid
- Ideal for quick spin
- **Includes:** both rotors and adapters for 0.5ml & 0.2ml • 1 YR WARRANTY

	M110941	M110942	M110943	M110944	M110945
Lid Color:	Blue	Clear	Neon Green	Purple	Red
Speed:	Fixed 6000 rpm / 2000 x g				
Capacity:	8 x 1.5 / 2.0 ml tubes; 16 x 0.2 ml PCR tubes; 2 x 0.2 ml PCR strips				
Dims:	4.5 x 5.9 x 4.5 in.				
Weight:	1.1 lbs.				
Accessories:					
M111116:	Replacement Adapters for 0.5 ml				
M111117:	Replacement Adapters for 0.2 ml				



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002