

## Mini Centrifuge, My Fuge™ 12



- Unlike traditional mini centrifuges, unit eliminates need to change rotors when switching between microtubes & PCR strips
- Runs (12) microtubes & (4) PCR strips simultaneously w/ the COMBI-Rotor™ (included)
- 5.5 Inches wide; near silent operation
- Starts/stops with open/close of lid
- Perfect for quick spin downs
- Quickly ramps up to 5500 rpm
- 1 YEAR WARRANTY

Cat. No.	Description
M110938	MyFuge <sup>12</sup>

### Adapters

M110939	Replacement Adapters for 0.5 ml
M110940	Replacement Adapters for 0.2 ml

	M110938
Speed:	Fixed 5500 rpm / 2000 x g
Capacity:	12 x 1.5/2.0 ml Tubes; 32 x 0.2 ml PCR Tubes; 4 x 0.2 ml PCR Strips
Dims. / Wgt.:	5.5 x 7.9 x 4.4 in. / 2.1 lbs.



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002