

Personal Mini-Centrifuge



- Powerful centrifugal force of 15100 g
- High speed, 12-Place Rotor is fully autoclavable and is made of high strength aluminum alloy
- Unique cooling design, minimizes rotor heating
- Over-speed detection; whisper quiet operation
- Brushless motor drive, maintenance free
- Pulse operation for quick spins
- Large LCD, auto door lock
- 2 YEAR WARRANTY

	M111319
Speed:	15000 rpm (200 - 15000rpm), 100 rpm Increments
Max. RCF:	15100 x g, 100 rpm x g Increments
Run Time:	30 secs. to 99 mins. ; Timer or Continuous
Dims./Wgt.:	Compact, 9.8 (W) x 9.6 (L) x 5.5 (H) in. / 13.2 lbs.



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002