

Portable pH Meter



Orion Star™

- AUTO-READ locks in the stable reading on screen & ready indicator alerts when readings are stable
- Up to 5 point pH calibration with auto recognition for NIST & DIN buffers
- Non-volatile memory holds up to 1000 datapoints with time and date stamp
- Portable, waterproof and protected to take anywhere; IP67-rated housing
- Complimentary data analysis software keeps meter software up-to-date
- Four AA batteries (included) for over 800 hours of operation or power adapter (sold separately)
- Large, backlit graphic LCD display shows pH/mV/RmV and temperature along with electrode status, time, date and calibration points
- Certifications: CE, TUV, FCC Class A limits
- 3 YEAR WARRANTY ON METER

M111357		
pH	Range:	-2.000 to 20.000
	Resol.:	0.1, 0.01, 0.001
	Relative Accuracy:	±0.002
Temp.	Range:	23° to 221°F (-5° to 105°C)
	Resol.:	0.1
	Relative Accuracy:	±0.1°C
mV	Range:	±2000.0
	Resol.:	0.1
	Relative Accuracy:	±0.2 mV or 0.05% (greater)
Relative Humidity:		5 to 85% Non-Condensing

pH Meter Kit includes: ROSS Gel-filled Epoxy pH/ATC Triode Electrode, pH 4, 7 and 10 buffers; Electrode Storage Solution; Rinse Solution Pouches; Hard Carry Case; HH Meter Armor with Stand and Electrode Holders



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002