

Compact Scale



- View in bright or dimly lit areas; 16mm LCD with amber backlight
- Counting mode with Automatic Counting Accuracy Improvement
- **Units of measure:** g, oz, lb, troy ozt, ct, mom, dwt., grains, newtons and teal
- 1 Second Stabilization; % Weighing
- **Features:** Over-under, specific gravity function and Stackable
- **Options:** RS-232 or USB, underhook, breeze break, carry case or a density determination kit
- AC power or battery (not included)
- 5 Year Warranty

	Cap. x Res.	Pan Size
M110972	120 g x 0.01 g	Ø4.3 in.
M110973	210 g x 0.01 g	Ø4.3 in.
M110974	310 g x 0.01 g	Ø4.3 in.
M110975	410 g x 0.01 g	Ø4.3 in.
M110976	610 g x 0.01 g	Ø4.3 in.
M110977	1500 g x 0.01 g	5 x 5.5 in.
M110978	2000 g x 0.01 g	5 x 5.5 in.
M110979	3000 g x 0.01 g	5 x 5.5 in.
M110980	4100 g x 0.01 g	5 x 5.5 in.
M110981	6100 g x 0.01 g	5 x 5.5 in.



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002