

Pinnacle Chek-Mite™



- 2 point calibration, Features manual or auto-calibration, auto-buffer recognition
- M102494 has auto temp. compensation
- M102492 has a manual hold button that freezes display
- M102493, M102494 and M102495 have Automatic End Point sensing (display holds the stable reading)
- M102495 has BNC connector for standard BNC electrode
- Battery (incl.) 500 hrs; waterproof case that float
- WARRANTY; 1 YR on Meter and 6 MONTH on Probe

	M102492
Range:	0 to 14 pH
Resolution:	0.01 pH
Weight:	6 oz.



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002