

# 3D Gyrotory Rockers, Lab Mart® Blotboy™



Mini BlotBoy



BlotBoy

- Optimized speed and tilt for gel blotting, washing and staining
- Tri-directional motion for thorough, gentle mixing
- Maintenance free, brushless motor; cold room and incubator safe
- Accepts stackable platforms and handles loads up to 2kg
- Bungee tie down points on platform, cords sold separately
- **Optional** stacking platform with flat mat, extra dimpled mat
- Units supplied with autoclavable flat mats
- 2 YEAR WARRANTY

Cat. No.	Description
M281000	Lab Mart® Mini BlotBoy™
M281010	Lab Mart® BlotBoy™

**Accessories**

M281001	Stacking Platform w/ Flat Mat for above
M281002	Dimpled Mat for above
M281003	Flat Mat 10.5 x 7.5 in.
M281004	Stacking Platform w/ Flat Mat for above
M281005	Dimpled Mat for above
M111118	Bungee Cord for both units

	M281000	M281010
<b>Speed:</b>	12RPM, Fixed	
<b>Motion Pitch:</b>	3D Motion / Fixed at 8°	
<b>Maximum Load:</b>	2 kg / 4.4 lbs. Centered and Balanced	
<b>Operating Temp. Range:</b>	4 to 65°C	
<b>Platform Size:</b>	10.5 x 7.5 in.	12 x 12 in.
<b>Dims. (WxDxH):</b>	10.5 x 7.5 x 8 in.	12 x 12 x 8 in.
<b>Weight:</b>	4.4 lbs.	



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002