

Mini-Stirrers, Magnetic



M280005



M280006

- Light, compact in ABS housing that resists most chemicals
- Comes with micro stir bar
- 2 YEAR WARRANTY

	M280005	M280006	M280007
Speed Range:	100 - 1000 RPM		
Environment:	0 to 50°, RH max 95%		
Dims.:	4.8 x 4.8 x 1.8 in.		
Cover:	ABS Plastic	Stainless Steel	ABS Plastic
Power:	110/220 VAC	110/220 VAC (Socket Supplied)	12V DC Car Battery
Weight:	1.4 lbs.	1.6 lbs.	1.3 lbs.

Cat. No.	Description
M280005	Ideal for Limited Bench Space in Lab
M280006	Creates Exothermic Reactions, AISI 316 M280007Field Unit, Low Input Voltage



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002